COVENTRY Cleanroom Wipes and Swabs

Coventry™ Sealed Foam Swabs

Economical versatile cleanroom swabs

The most cost-effective swabs for general-purpose applications. 100 ppi reticulated medical-grade polyurethane foam provides excellent particle entrapment. The open-cell structure gives the best solvent holding capacity of any style of swab. These swabs are free of adhesives or binders that could contaminate your process. Fiberless construction does not generate particles.



| Product | Pkg Count | Swab Length | Handle Width | Head Length | Head Width | Handle Material | Head Material |
|-------------------------|--------------------|----------------|---------------------|--------------------|----------------------------------|------------------------|----------------------|
| 41050 | 500 swabs | 2.8" | 0.10" | 0.40" | 0.13" x 0.14" | polypropylene | polyurethane |
| | | 7.2 cm | 2.6 mm | 10.2 mm | 3.3 mm x 3.5 mm | | foam |
| | | | | | | | |
| 41050I | 500 swabs | 2.8" 7.2 cm | 0.10" 2.6 mm | 0.40" 10.2 mm | 0.13" x 0.14" 3.3 mm x 3.5 mm | polypropylene | polyurethane foam |
| AQUA | | | 2.0 11111 | 10.2 11111 | 0.0 mm x 0.0 mm | | ioam |
| PRIME | or extra cleanline | ess: | | | | | |
| | | | | | | | |
| 43170 | 500 swabs | 6.0" | 0.09" | 0.80" | 0.33" x 0.35" | wood | polyurethane |
| | | 15.2 cm | 2.2 mm | 20.3 mm | 8.4 mm x 8.8 mm | | foam over cotton |
| | | | | | | | |
| | | | | | | | |
| | | | | | | _ | |
| 14070 | 500 swabs | 3.2" | 0.10" | 0.39" | 0.11" x 0.13" | polypropylene | polyurethane |
| 44070 | 500 swabs | 3.2" 8.0 cm | 0.10" 2.4 mm | 0.39" 9.9 mm | 0.11" x 0.13" 2.7 mm x 3.3 mm | polypropylene | polyurethane foam |
| 44070 | 500 swabs | | | | | polypropylene | |
| | 500 swabs | 8.0 cm 2.7" | 2.4 mm 0.09" | 9.9 mm 0.41" | 2.7 mm x 3.3 mm 0.10" x 0.12" | polypropylene | foam polyurethane |
| | | 8.0 cm | 2.4 mm | 9.9 mm | 2.7 mm x 3.3 mm | | foam |
| | | 8.0 cm 2.7" | 2.4 mm 0.09" | 9.9 mm 0.41" | 2.7 mm x 3.3 mm 0.10" x 0.12" | | foam polyurethane |
| 44070 48040 49280 | | 8.0 cm 2.7" | 2.4 mm 0.09" | 9.9 mm 0.41" | 2.7 mm x 3.3 mm 0.10" x 0.12" | | foam polyurethane |



PRIME For Extra Cleanliness

When cleaned with the Aqua Prime process, Coventry Sealed Foam Swabs minimize background contamination making them suitable for both critical cleaning and cleaning validation procedures. Aqua-Prime cleaning process can be added to any swab with special order.